

# **RABBIT™ RESTRAINER**

- **Safety for the animal**
- **Practicality of use**
- **Flexibility of use**



**HOUSING**

 **TECNIPLAST**  
*innovation through passion<sup>SM</sup>*

Easy to use for pyrogen testing, injection or auricular bleeding, while providing maximum security and protection. Only for rabbits weighing from 3 to 4,5 Kg.

### 1/Safety for the animal

Specially designed top plate **protects the animal's spine**, inverted U-shaped grid provides **paw comfort** and neck collar to protect the animal during procedure.

### 2/Practicality of use

Stainless steel (S/S) base insert (267 x 108 mm) provides firm footing while allowing faeces and urine to drop into the tray below. Bottom plate is slotted at the rear for pyrogen testing. All parts are **autoclavable and easy to clean**.

### 3/Flexibility of use

Rack for restrainer is also available with the second tier located 1460mm from floor level.




## HOUSING

### Rabbit Restrainer

Safety for the animal / Practicality of use

**1&2**



Flexibility of Use

**3**



#### TECHNICAL SHEET

	Dimensions (L x W x H) mm / in
Restrainer* (PC)	400 x 200 x 200 / 15.75 x 7.87 x 7.87
Tray (PC)	298 X 127 X 20 / 15.75 x 7.87 x 7.87

\*Waste tray not included.

PC = Polycarbonate.

No. of Cages per Rack	Dimensions (L x W x H) mm / in
6	1475 x 500 x 750 / 58.07 x 19.68 x 29.53
9	2190 x 500 x 750 / 86.22 x 19.68 x 29.53
12	1475 x 500 x 1250 / 58.07 x 19.68 x 49.21



ISO-9001  
Cert. n° 0875



ISO-14001  
Cert. n° 0256A/0



HEADQUARTERS  
ITALY | Tel. +39 0332 809711 • www.tecniplast.it • E-mail: tecnicom@tecniplast.it

SISTER COMPANIES

AUSTRALIA/NEW ZEALAND | Tel. + 61 2 8845 6500 • www.tecniplast.it • E-mail: info@tecniplast.com.au

CHINA | Tel. + 86 21 50810920 • www.tecniplast.cn • E-mail: tecniplastchina@tecniplast.it

FRANCE | Tél. 04 72 52 94 41 • www.tecniplast.fr • E-mail: info@tecniplast.fr

GERMANY | Tel. 08805 921320 • www.tecniplast.de • E-mail: info@tecniplast.de

UNITED KINGDOM | Tel. 0845 0504556 • www.tecniplast.it • E-mail: info@tecniplastuk.com

JAPAN | Tel. + 03 5770 5375 • www.tecniplast.it • E-mail: info@tecniplastjapan.co.jp

USA/CANADA | Toll Free: 877.669.2243 • www.tecniplastusa.com • E-mail: info@tecniplastusa.com

To find your local distributor, please visit [www.tecniplast.it](http://www.tecniplast.it)

Tecniplast reserves the right to change or modify product and or specifications without notice or obligation. Rabbit Restrainer is a Tecniplast trademark.